





Meet the New Kid on the Block

- Smooth Roll-Through
- Dynamic Toe-Off
- Most Durable Low-Profile Foot Available



TECHNICAL SPECS

Item N°.	KID
Weight (26 Cat 4, with foot shell and spectra sock)	505g
Build Height*	2" / 51mm
Size	19 - 21cm
Heel Height	3/8" / 10mm

*Build height includes spectra sock, foot shell and 10mm heel height.


Black Spectra Sock and Heel Wedge included in each Foot Module delivery.
Foot Shell must be ordered separately.

WEIGHT lbs./ kg.	IMAPCT LEVEL FOR SIZES 19-21cm			
	LOW	MODERATE	HIGH	HIGH IMPACT
0-30lbs / 0-14kg	1P	1P	1P	2P
31-45lbs / 15-20kg	2P	2P	2P	3P
46-60lbs / 21-27kg	3P	3P	3P	4P
61-75lbs / 28-34kg	4P	4P	4P	5P
76-90lbs / 35-41kg	5P	5P	5P	6P
91-105lbs / 42-48kg	6P	6P	6P	7P
106-120lbs / 49-54kg	7P	7P	7P	8P
131-135lbs / 55-61kg	8P	8P	8P	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

Complete Prosthetic Foot Part Number

Foot Module Part Number



*Category selections vary with sizes. Please refer to category selections available for each product.

Foot Shell Part Number

