





Virtually Indestructible

- Highly-dynamics
- High-energy return
- Smooth roll-over
- All-terrain and waterproof (fresh or salted water)



TECHNICAL SPECS

Item N°.	HiP	
Weight (26 Cat 4, with foot shell and spectra sock)	920g	
Build Height*	22 - 24cm	6.5" / 165mm
	25 - 27cm	6.875" / 175mm
	28 - 29cm	7.125" / 181mm
Heel Height	3/8" / 10mm	

*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell and 10mm heel height.

WEIGHT lbs./ kg.	IMPACT LEVEL			
	LOW	MODERATE	HIGH	HIGH IMPACT
0-105lbs / 0-48kg	1	1	1	2
106-140lbs / 49-64kg	2	2	2	3
141-175lbs / 65-79kg	3	3	3	4
176-210lbs / 80-95kg	4	4	4	5
211-245lbs / 96-111kg	5	5	5	6
246-280lbs / 112-127kg	6	6	6	7
281-315lbs / 128-143kg	7	7	7	8
316-350lbs / 144-159kg	8	8	8	9
351-365lbs / 160-166kg	9	9	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.




L5981

Also available in EVAQ8 and H2O Models

Black Spectra Sock and Heel Wedge included in each Foot Module delivery.
Foot Shell must be ordered separately.

Complete Prosthetic Foot Part Number

Foot Module Part Number



*Category selections vary with sizes. Please refer to category selections available for each product.

Foot Shell Part Number

