





Unbridled Capabilities

- Cosmesis-ready to meet your aesthetic preferences
- Lower overall product weight
- Smooth roll through with unparalleled energy return
- Biometric ankle motion simulates natural gait, increasing patient comfort
- Extreme reliability for any use in any environment



Also available in EVAQ8 and H2O Models



TECHNICAL SPECS


Item N°	RAM	
Weight (26 Cat 4, with foot shell and spectra sock)	850g	
Build Height*	22 - 24cm	6.625" / 168mm
	25 - 27cm	7.125" / 181mm
	28 - 29cm	7.25" / 184mm
Heel Height	3/8" / 10mm	

*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell and 10mm heel height.

Black Spectra Sock and Heel Wedge included in each Foot Module delivery. Foot Shell must be ordered separately.


Complete Prosthetic Foot Part Number

Foot Module Part Number



*Category selections vary with sizes. Please refer to category selections available for each product.

Foot Shell Part Number



WEIGHT lbs./kg.	IMPACT LEVEL			
	LOW	MODERATE	HIGH	HIGH IMPACT
0-105lbs / 0-48kg	1	1	1	2
106-140lbs / 49-64kg	2	2	2	3
141-175lbs / 65-79kg	3	3	3	4
176-210lbs / 80-95kg	4	4	4	5
211-245lbs / 96-111kg	5	5	5	6
246-280lbs / 112-127kg	6	6	6	7
281-315lbs / 128-143kg	7	7	7	8
316-350lbs / 144-159kg	8	8	8	9
351-365lbs / 160-166kg	9	9	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.