





Innovation You Can Feel

- Low-profile build
- Smooth roll-over
- Cosmesis-friendly
- Waterproof (fresh or salt water)
- Superior ground compliance
- High-energy return
- All-terrain
- Reduced socket issues




TECHNICAL SPECS

Item N°	RAMLP	
Weight (26 Cat 4, with foot shell and spectra sock)	805g	
Build Height*	22 - 24cm	4.5" / 114mm
	25 - 27cm	4.75" / 121mm
	28 - 29cm	4.75" / 121mm
Heel Height	3/8" / 10mm	

*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell and 10mm heel height.

WEIGHT lbs./ kg.	IMAPCT LEVEL			
	LOW	MODERATE	High	High Impact
0-105lbs / 0-48kg	1	1	2	3
106-140lbs / 49-64kg	2	2	3	4
141-175lbs / 65-79kg	3	3	4	5
176-210lbs / 80-95kg	4	4	5	6
211-245lbs / 96-111kg	5	5	6	7
246-280lbs / 112-127kg	6	6	7	8
281-315lbs / 128-143kg	7	7	8	9
316-350lbs / 144-159kg	8	8	9	-
351-365lbs / 160-166kg	9	9	-	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.



Also available in EVAQ8 and H2O Models

Black Spectra Sock and Heel Wedge included in each Foot Module delivery.
Foot Shell must be ordered separately.

Complete Prosthetic Foot Part Number

Foot Module Part Number



*Category selections vary with sizes. Please refer to category selections available for each product.

Foot Shell Part Number

