



1M SERIES




KNEES

With its convenient lock, the 1M Series ensures patient safety while walking.




- Joint friction adjustment.
- Lock with adjustable stop.
- Socket rotation adjustment: $\pm 15^\circ$.

1M01 / 1M10



1M01-P6 / 1M10-P6



TECHNICAL FEATURES

| Item N°. | Max. Patient's Weight | Weight | TH | BH | H | Material | Release lever | Flexion | Distal Adapter |
|----------|-----------------------|--------|--------|-------|-------|----------------|---------------|---------|-----------------|
| 1M01 | 220 lbs/ 100 kg | 345 g | 99 mm | 58 mm | 13 mm | Alu + Steel | No | 130 ° | Tube Ø 30 mm |
| 1M10 | | 355 g | | | | | Yes | | |
| 1M01-P6 | 275 lbs/ 125 kg | 375 g | 102 mm | 59 mm | 14 mm | | No | | Tube Ø 34 mm |
| 1M10-P6 | | | | | | | Yes | | |

ORDER ONE OF THE FOLLOWING TO COMPLETE YOUR ORDER



1K40

1K160-P6

1K160

1K163

ASK A QUESTION

ORDER