

RUSH FOOT® EVAQ8 COLLECTION



Simplistic, Effective, Efficient


- Simple, integrated design
- Stabilizes volume, increases linkage, and heightens proprioception through elevated vacuum suspension
- All-terrain, fully-submersible
- Requires minimal maintenance
- The RUSH Foot® EVAQ8 Collection is available in the RAMPAGE, HiPro, ROGUE 2, and RAMPAGE LP models




L5987, L5984, L5781

Complete Prosthetic Foot Part Number

Foot Module Part Number



Foot Shell Part Number





*Category selections vary with sizes. Please refer to category selections available for each product.

Rebuild KIT (EVRB)

- 90 Degree Hose Barb
- 2' Tygon Vacuum Hose
- Hose Barb Filter
- 1/2" Wide Hook and Loop

1X998
To avoid spoiling
the pyramid
during trials



TECHNICAL SPECS

Item N°.		RAMPAGE: EVRAM	HiPro: EVQH	RAMPAGE LP: EVQRAML P	ROGUE: EVQR2
Weight (26 Cat 4, with foot shell and spectra sock)		860g	930g	815g	1048g
Build Height*	22 - 24cm	6.625" / 168mm	6.5" / 165mm	4.5" / 114mm	6.125" / 155mm
	25 - 27cm	7.125" / 181mm	6.875" / 175mm	4.75" / 121mm	6.625" / 168mm
	28 - 29cm	7.25" / 184mm	7.125" / 181mm	4.75" / 121mm	6.750" / 171mm
Heel Height		3/8" / 10mm			

*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell and 10mm heel height.

Black Spectra Sock and Heel Wedge included in each Foot Module delivery.
Foot Shell must be ordered separately.