




## Extraordinary Adventures Begin with the ROVER

- Smooth Roll-Through
- Dynamic Toe-Off
- Most Durable Low-Profile Foot Available



L5981



### TECHNICAL SPECS

Item N°.	ROV
Weight (26 Cat 4, with foot shell and spectra sock)	950g
Build Height*	22 - 24cm 2.125" / 54mm
	25 - 27cm 2.375" / 60mm
	28 - 30cm 2.5" / 64mm
Heel Height	3/8" / 10mm

\*Build height based on size 23cm, 26cm, or 29cm Cat 4 foot module with spectra sock, foot shell and 10mm heel height.


Black Spectra Sock and Heel Wedge included in each Foot Module delivery.  
Foot Shell must be ordered separately.

WEIGHT lbs./ kg.	IMAPCT LEVEL			
	LOW	MODERATE	High	High Impact
0-105lbs / 0-48kg	1	1	1	2
106-140lbs / 49-64kg	2	2	2	3
141-175lbs / 65-79kg	3	3	3	4
176-210lbs / 80-95kg	4	4	4	5
211-245lbs / 96-111kg	5	5	5	6
246-280lbs / 112-127kg	6	6	6	7
281-315lbs / 128-143kg	7	7	7	8
316-350lbs / 144-159kg	8	8	8	9
351-365lbs / 160-166kg	9	9	9	-

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability.

### Complete Prosthetic Foot Part Number

#### Foot Module Part Number



#### Foot Shell Part Number

