

FREEDOM

Pacifica®

FREEDOM

Pacifica® LP




Because Height Matters

- Water-resistant
- Choose between the **Pacifica® (FS2)** and **Pacifica® LP (FS4)**, each perfect for amputees with long residual limbs or those needing more build space
- LP design requires as little as 2" of clearance
- Lightweight and durable
- Rated to 365lbs, 36-month warranty



Sandal-toe available for sizes 22-28cm.



TECHNICAL SPECS

Item N°.	FS2: Pacifica® (sizes: 24 - 31cm)	FS4: Pacifica® LP (sizes: 22 - 31cm)
Weight (27 Cat 4, module with foot shell)	545g	537g
Build Height*	22 - 25cm	2.5" / 65mm
	26 - 28cm	2.6" / 66mm
	29 - 31cm	2.7" / 70mm
Heel Height	3/8" / 10mm	3/8" / 10mm


*Build height based on size 23cm, 26cm, 29cm Cat 4 module with foot shell and 10mm heel height. Sandal-toe available for sizes 22-28cm.

Black Spectra Sock and Heel Bumper Kit included in each Foot Module delivery.
Foot Shell must be ordered separately.

WEIGHT lbs./ kg.	IMAPCT LEVEL		
	LOW	MODERATE	HIGH
100-115lbs / 44-52kg	1	1	2
116-130lbs / 53-59kg	1	2	3
131-150lbs / 60-68kg	2	3	4
151-170lbs / 69-77kg	3	4	5
171-195lbs / 78-88kg	4	5	6
196-220lbs / 89-100kg	5	6	7
221-255lbs / 101-116kg	6	7	8
256-285lbs / 117-130kg	7	8	9
286-325lbs / 131-147kg	8	9	-
326-365lbs / 148-166kg	9	-	-

Complete Prosthetic Foot Part Number

Foot Module Part Number



Foot Shell Part Number

